




<https://www.fairmat-nfdi.eu/fairmat/consortium>



NOVEL MATERIALS DISCOVERY

<https://nomad-lab.eu/>

Ellipsometry Example and Tutorial in NOMAD



Create and account or login using an existing account

NOMAD – Manage and Publish Materials Data

This is the *graphical user interface* (GUI) of NOMAD. It allows you to **search, access, and download all NOMAD data** in its *raw files* and *processed data* form. You can **upload and manage your own raw materials science data**. You can access all published data without an account. If you want to provide your own data, please login or register for an account.

You can learn more about NOMAD on its [homepage](#), our [documentation](#). There is also an [FAQ](#) and the more detailed [uploader documentation](#).

Interactive Search

NOMAD extracts **rich metadata** from uploaded raw-data. [Explore NOMAD's data](#) by creating complex queries from interactive data visualizations of key properties, including the simulated composition/system, used method, upload metadata, as well as material classifications and available quantities. Or use the **Optimade** filter language to add arbitrarily nested queries.

A common data format

NOMAD provides data in *processed* and *normalized* form in a machine processable and common hierarchical format. This *processed data*, i.e. the **NOMAD Archive**, is organized into nested sections of quantities with well defined units, data types, shapes, and descriptions. These definitions are called the **NOMAD Metainfo** and they can be [browsed here](#).

What is shown on these slides is a developer's version of NOMAD, but the example will be integrated into NOMAD soon

You can create an upload and upload files through this browser-based interface:

[CREATE A NEW UPLOAD](#) or [ADD EXAMPLE UPLOADS](#)

Or, you can create an upload by sending a file-archive via shell command:

```
curl -X POST 'http://localhost:8000/fairdi/nomad/latest/api/v1/uploads?token=A2wS54YUQzKqSxw13V23Tw.tLL7tjTKTi4jQqPSFyzXXlf9ZBc' -T <local_file>
```

In PUBLISH → UPLOADS either a new upload can be created or an existing example can be uploaded

You can create an upload and upload files through this browser... CREATE A NEW UPLOAD or ADD EXAMPLE UPLOADS Or, you can create an upload by sending a file-archive via shell... curl -X POST 'http://localhost:8000/fairdi/nomad'

Select a sample Upload

Electronic structure code input and output files ADD

This upload demonstrate the basic use of NOMAD's *parsers*. For many *electronic structure codes* (VASP, etc.), NOMAD provides parsers. You simply upload the *input and output files* of your simulations and NOMAD parsers are extracting all necessary metadata to produce a **FAIR** dataset.

Electronic Lab Notebook ADD

This example contains a custom NOMAD *schema* to create an **Electronic Lab Notebook (ELN)** and a few example *data* entries that use this schema. The schema demonstrates the basic concepts behind a NOMAD ELN and can be a good **starting point** to create you own schemas that model **FAIR data** acquired in your lab.

Tabular Data ADD

This upload demonstrates the used of tabular data. In this example we use an *xlsx* file in combination with a custom schema. The schema describes what the columns in the excel file mean and NOMAD can parse everything accordingly to produce a **FAIR** dataset.

Ellipsometry ADD




This example presents the capabilities of the NOMAD platform to

CANCEL



OVERVIEW

FILES

 **Ellipsometry** 
upload id: oKltn96XQSieCGNTwkHn7Q 



Each upload has its own upload id

README.md

Introduction

This example presents the capabilities of the NOMAD platform to store and standardize ellipsometry data. It shows the generation of a NeXus file according to the [NXellipsometry](#) application definition and a successive analysis of an example data set (measured values of the ellipsometric angles Psi and Delta of SiO2 on Si).

Viewing uploaded data

Below, you find an overview of your uploaded data. Click on the `> /` button to get a list of your data or select **FILES** from the top menu of this upload. You may add your own files to the upload or experiment with the pre-existing electronic lab notebook (ELN) example. The ELN follows the general structure of NOMAD ELN templates and you may refer to the [documentation](#) or a [YouTube tutorial](#) (~1h) for further information. When the ELN is saved a NeXus file will be generated from the provided example data. You may also view your supplied or generated NeXus files here with the H5Web viewer. To do so open the **FILES** tab and just select a `.nxs` file.



OVERVIEW

FILES

DATA

LOGS

Metadata

type Ellips name ellips.data.archive.json comment no comment references authors Carola Emminger datasets no datasets mainfile ellips.data.archive.json entry id O3r5hUaa8xjPvPt8im_kJ3o2LsL_ material id unavailable upload id oKltn96XQSieCGNTwkHn7Q

Ellips

Reader ellips NxdI NXellipsometry.nxdI input_files test-data.dat input_files eln_data.yaml input_files Output SiO2onSi.ellips.nxs Filename test-data.dat

Electronic Lab Notebook (ELN)

Metadata according to Nxellipsometry
(https://manual.nexusformat.org/classes/contributed_definitions/NXellipsometry.html)

Light source
arc lamp

Focussing probes

Angular spread
0.2 Unit rad^2

Ellipsometry type
dual compensator

Calibration status
no calibration

Stage type
manual stage

Detector type
CCD spectrometer

Count time
1 Unit s

Integration time
0 Unit s

Rotating element
compensator (source side)

Atom types
Si, O

Electronic Lab Notebook (ELN)

All files contained in the example upload can be found in FILES



- Entry files**
- Ellipsometry workflow example.ipynb
 - README.md
 - SiO2onSi.ellips.nxs
 - Si_Aspnes.mat
 - ellips.data.archive.json nomad
 - ellips.scheme.archive.yml nomad
 - eln_data.yaml**
 - test-data.dat

eln_data.yaml

file

file size
1.7 kB (1688 bytes)

preview

```
acquisition_program/\@url: https://www.jawoollam.com/ellipsometry-software/completeease
acquisition_program/program: CompleteEASE
acquisition_program/version: '6.37'
address: Zum Großen Windkanal 2, 12489 Berlin, Germany
affiliation: Humboldt-Universität zu Berlin
angle_of_incidence_unit: degrees
angular_spread:
  unit: rad ** 2
  value: 0.2
atom_types: Si, O
blocks:
- type
- angle_of_incidence
calibration_status: no calibration
colnames:
- type
- wavelength
- angle_of_incidence
- psi
- delta
- err.psi
- err.delta
column_names:
- psi
- delta
company: J. A. Woollam Co.
count_time:
  unit: s
  value: 1
data_identifier: 1
data type: psi/delta
```

When saving the ELN, a YAML file is created that lists all the metadata defined in the ELN



OVERVIEW

FILES

DATA

LOGS

Entry files

- folder
- Ellipsometry workflow example.ipynb
- README.md
- SiO2onSi.ellips.nxs
- Si_Aspnes.mat
- ellips.data.archive.json nomad
- ellips.scheme.archive.yml nomad
- eln_data.yaml
- test-data.dat**

test-data.dat

file

file size

482 kB (482497 bytes)

preview

2nm SiO2 on Si on RC2

VASEmethod[EllipsometerType=4 , CompleteEASE=6.37, AcqTime=15.000, ZoneAve=1, Acq. Parameters=DE

Angstroms

E	1930.000000	50.000000	40.014217	142.127655	0.008585	0.034774
E	1940.000000	50.000000	40.026409	142.307449	0.007827	0.031627
E	1950.000000	50.000000	40.031204	142.483231	0.007221	0.029106
E	1960.000000	50.000000	40.045902	142.679901	0.006813	0.027399
E	1970.000000	50.000000	40.065468	142.866608	0.006478	0.025995
E	1980.000000	50.000000	40.077728	143.102463	0.006224	0.024925
E	1990.000000	50.000000	40.080673	143.409317	0.006016	0.024045
E	2000.000000	50.000000	40.082367	143.746017	0.005829	0.023258
E	2010.000000	50.000000	40.072388	143.990952	0.005649	0.022496
E	2020.000000	50.000000	40.073170	144.279099	0.005475	0.021767
E	2030.000000	50.000000	40.069405	144.611435	0.005311	0.021084
E	2040.000000	50.000000	40.066757	144.883820	0.005150	0.020417
E	2050.000000	50.000000	40.061028	145.150055	0.004993	0.019760
E	2060.000000	50.000000	40.055706	145.396118	0.004843	0.019128
E	2070.000000	50.000000	40.036533	145.661514	0.004703	0.018531
E	2080.000000	50.000000	40.017166	145.896454	0.004572	0.017967
E	2090.000000	50.000000	40.025864	146.114670	0.004443	0.017414
E	2100.000000	50.000000	40.028156	146.370560	0.004322	0.016892
E	2110.000000	50.000000	40.012478	146.607407	0.004227	0.016464
E	2120.000000	50.000000	40.007420	146.799438	0.004136	0.016055
E	2130.000000	50.000000	40.005928	147.055832	0.004045	0.015653
E	2140.000000	50.000000	40.001190	147.294281	0.003968	0.015308
E	2150.000000	50.000000	39.999153	147.531891	0.003896	0.014983
E	2160.000000	50.000000	39.998287	147.757156	0.003830	0.014690
E	2170.000000	50.000000	39.997860	147.958969	0.003769	0.014416
E	2180.000000	50.000000	40.001450	148.183533	0.003709	0.014154
E	2190.000000	50.000000	39.995445	148.426010	0.003648	0.013888
E	2200.000000	50.000000	39.998291	148.655228	0.003589	0.013635

The file containing the measured data



OVERVIEW

FILES

DATA

LOGS

- Entry files
- Ellipsometry workflow example.ipynb
 - README.md
 - SiO2onSi.ellips.nxs
 - Si_Aspnes.mat
 - ellips.data.archive.json
 - ellips.scheme.archive.yml
 - eln_data.yaml
 - test-data.dat

Ellipsometry workflow example.ipynb

file

file size
164 kB (163834 bytes)

TOOLS

jupyter: Basic jupyter run with an empty notebook or on given notebook file.
LAUNCH

nexustools: Includes multiple NeXus tools for visualization and analysis.
LAUNCH

ellips: An example for analyzing ellipsometric data.
LAUNCH

mpes: An example for analyzing mpes data.

Launch "ellips" in jupyter lab to view the tutorial "Ellipsometry workflow example"

Filter files by name

/ uploads / oKltn96XQSieCGNTwkHn7Q /

Name	Last Modified
ellips.data.archive.json	40 minutes ago
ellips.scheme.archive.yml	40 minutes ago
Ellipsometry workflow example.ipynb	9 minutes ago
eln_data.yaml	40 minutes ago
README.md	40 minutes ago
Si_Aspnes.mat	40 minutes ago
SiO2onSi.ellips.nxs	40 minutes ago
test-data.dat	40 minutes ago

The tutorial explains how to create and inspect the NeXus file using h5web

Launcher Ellipsometry workflow exam

Code Python 3 (ipykernel)

Ellipsometry workflow example

In this notebook, an ellipsometry data set of 2 nm SiO₂ on Si is analyzed using the analysis tool [pyElli](#)

1. Create NeXus file from measurement data

The metadata of the experiment are listed in a YAML file (**eln-data.yaml**, which is automatically created when saving the metadata entered into the electronic lab notebook (ELN) within NOMAD) according to the application definition [NXellipsometry](#). The name of the data file (here **test-data.dat**) needs to be specified in the ELN and, hence, is defined as an entry 'filename' in the YAML file. Using the **ellips** reader and the application definition in NXDL format, a NeXus file (**SiO2onSi.ellips.nxs**) is created. Both the data and metadata files must be stored in this repository.

Note: When creating or modifying the YAML file without using the ELN, make sure that all required fields are provided; recommended and optional fields may be provided if known and meaningful.

```
[1]: from nexusparser.tools.dataconverter.convert import convert
```

```
[ ]: convert(input_file=["eln_data.yaml"],
            reader='ellips',
            nxdl='NXellipsometry',
            output='SiO2onSi.ellips.nxs')
```

2. Inspect the NeXus file with h5web

```
[ ]: from jupyterlab_h5web import H5Web
```

File Edit View Run Kernel Tabs Settings Help

Filter files by name

/ uploads / oKltn96XQSieCGNTwkHn7Q /

Name	Last Modified
ellips.data.archive.json	41 minutes ago
ellips.scheme.archive.yml	41 minutes ago
Ellipsometry workflow example.ipynb	10 minutes ago
eln_data.yaml	41 minutes ago
README.md	41 minutes ago
Si_Aspnes.mat	41 minutes ago
SiO2onSi.ellips.nxs	seconds ago
test-...dat	41 minutes ago

Name: SiO2onSi.ellips.nxs
 Size: 90.2 KB
 Path: uploads/oKltn96XQSieCGNTwkHn7Q
 Created: 2022-09-21 17:30:04
 Modified: 2022-09-21 17:30:04
 Writable: true

NeXus file: contains all the metadata and the measured data (+plots)

Launcher Ellipsometry workflow exam

Python 3 (ipykernel)

2. Inspect the NeXus file with h5web

```
[3]: from jupyterlab_h5web import H5Web
```

```
[4]: H5Web('SiO2onSi.ellips.nxs')
```

```
[4]:
```

SiO2onSi.ellips.nxs

entry NX

Display Inspect Feedback

NX Spectrum Linear Linear Auto-scale

psi_50deg (degrees)

psi_50deg (degrees)

wavelength (angstrom)

This is the end of the general template. Continue to fill the notebook based on **your own** post-processing of the *.nxs file.

Filter files by name

/ uploads / oKltn96XQSieCGNTwkHn7Q /

Name	Last Modified
ellips.data.archiv...	42 minutes ago
ellips.scheme.ar...	42 minutes ago
Ellipsometry wor...	seconds ago
eln_data.yaml	42 minutes ago
README.md	42 minutes ago
Si_Aspnes.mat	42 minutes ago
SiO2onSi.ellips....	a minute ago
test-data.dat	42 minutes ago

- entry NX
 - acquisition_program NX
 - program
 - version
 - definition
 - experiment_description**
 - experiment_identifier
 - instrument
 - operator
 - plot NX
 - sample
 - start_time

entry > experiment_description

Scalar

RC2 scan on 2nm SiO2 on Si in air

Display Inspect Feedback

These are the metadata as defined in the ELN shown above in accordance with NXellipsometry

Filter files by name

/ uploads / oKltn96XQSieCGNTwkHn7Q /

Name	Last Modified
ellips.data.archiv...	42 minutes ago
ellips.scheme.ar...	42 minutes ago
Ellipsometry wor...	a minute ago
eln_data.yaml	42 minutes ago
README.md	42 minutes ago
Si_Aspnes.mat	42 minutes ago
SiO2onSi.ellips....	2 minutes ago
test-data.dat	42 minutes ago

Launcher

Ellipsometry workflow exam... SiO2onSi.ellips.nxs

SiO2onSi.ellips.nxs

- entry NX
 - acquisition_program NX
 - program
 - version
 - definition
 - experiment_description
 - experiment_identifier
 - instrument
 - angle_of_incidence
 - angular_spread
 - calibration_status
 - company
 - detector
 - detector_type**
 - rotating_element
 - ellipsometry_type
 - firmware

entry > instrument > detector > **detector_type** Display Inspect Feedback

<> Scalar

CCD spectrometer

NOMAD JupyterLab

File Edit View

Filter files by

/ uploads / oKltn96XQSie

Name	Time
ellips.data.archiv...	an hour ago
ellips.scheme.ar...	an hour ago
Ellipsometry wor...	6 minutes ago
eln_data.yaml	an hour ago
README.md	an hour ago
Si_Aspnes.mat	an hour ago
SiO2onSi.ellips....	6 minutes ago
test-data.dat	an hour ago

The data (in this example, these are psi and Delta at three different angles of incidence) are contained in measured_data and can be viewed in form of a matrix, a line (shown here) or a heatmap

- > instrument
- > operator
- > plot NX
- ▼ sample
 - ⊞ atom_types
 - ⊞ column_names
 - ⊞ data_identifier
 - ⊞ data_type
 - ▼ environment_conditions
 - ⊞ medium
 - ⊞ layer_structure
 - ⊞ **measured_data**
 - ⊞ sample_history
 - ⊞ sample_name
 - ⊞ start_time

entry > sample > measured_data

Matrix Line Heatmap

n	1	1	3	2	1088
x	D0	D1	D2	D3	D4
D0	0:0	0:0	0:2	0:1	

measured_data

SiO2onSi.ellips.nxs

Questions or Comments?

If you have questions or comments regarding the example, please feel free to contact us:

Carola Emminger: emminger.carola@physik.hu-berlin.de

Florian Dobener: florian.dobener@physik.hu-berlin.de

